

**SAVE
20%**
on your medical
equipment
calibration!

Send us a copy of your previous invoice for your medical equipment calibration and/or PAT testing and we promise to beat the cost by 20% when you sign up for a visit from Calibrate UK. This offer applies to all new Calibrate UK customers.

1) Take a photo/scan of your invoice for your latest medical equipment calibration/PAT testing visit.

2) Email the invoice to info@calibrate-uk.co.uk and quote SAVE20% in the subject box.

3) A member of our team will get in touch to put a date in the diary!

Why switch to Calibrate UK?

- ✓ Fair prices from the start.
- ✓ Combine medical equipment calibration, PAT testing of non-medical items, LOLER inspections, Vaccine fridge servicing all in one visit.
- ✓ Professional and friendly technicians.
- ✓ Expert repairs for a range of manufacturers.

Quote form on reverse ⇨

For current customers: if you are already a Calibrate UK customer please fill out the section below and then pass this flyer to another practice manager.. For any successful referrals we will send you a FREE TEMPERATURE DATALOGGER for detailed monitoring of your Vaccine Fridge.

Surgery name:

Postcode:

Contact name:



QUOTATION FORM

*****FAO: THE PRACTICE MANAGER Re. Medical Equipment Testing & Calibration*****
Tel: 0845 643 5253 Fax: 0844 358 3365 Email: info@calibrate-uk.co.uk

For a no obligation quotation please fill out and fax return the form below:

Customer Details:
Practice name:
Contact name:
Practice Address:
Practice postcode:
Telephone number:
Fax number:
Email address:
Month Calibration due:

Medical Equipment details:

Item	No. of items		No. of items
Thermometers		Nebulisers	
Vaccine Refridgerators		Spirometers	
Height Measures		Medical Scales	
Hand Held Ultrasounds		Pulse Oximeters	
Auroscopes/ Oscopes		Blood Pressure Monitors	
Defibrillators		Syringe drivers/ feed pumps	
ECG Machines		Hearing loops	
Cauteries/ Hyfrecators		Blood Glucose Monitors	
Audiometers		Suction Pumps/ Aspirators	
Smokerlyzers (CO meters)		Other:	
Ear syringes		Electrical items for PAT testing	

FAX BACK TO 0844 358 3365
or scan to info@calibrate-uk.co.uk